





# Rubber Dam Clamp Forceps



## CLAMPS ( PLIERS )





# Rubber Dam Clamps

Anterior



**2510 - 9**  
Winged  
Labial Clamp



**2511 - 9T**  
Winged  
Edged Labial Clamp



**2512 - 212**  
Winged  
Cervical Clamp



**2513 - 210**  
Winged  
Labial Clamp

Premolars



**2514 - 00**  
Winged



**2515 - 0**  
Winged



**2516 - W0**  
Wingless

For Small Teeth, Elongated or With Narrow Collars



**2517 - 1**  
Winged



**2518 - 1T**  
Winged



**2519 - W1**  
Wingless

For General Use With Adjustable Jaws According  
To Gingival Line Shapes



**2520 - 2A**  
Winged



**2521 - 2AT**  
Edged, Winged



**2522 - W2A**  
Wingless

For Premolars, With Elongated Collar



**2523 - 2**  
Winged



**2524 - 2T**  
Edged, Winged



**2525 - W2**  
Wingless

For Premolars, With Elongated Collar



**2526 - 29**  
Wingless



**2527 - 60**  
Wingless



# Rubber Dam Clamps



**2528 - 3**  
Winged

**2529 - W3**  
Wingless

**2530 - 4**  
Winged

**2531 - W4**  
Wingless

**For Small Lower Molars**

**For Small Lower Molars With Round clamps**



**2532 - 7**  
Winged

**2533 - W7**  
Wingless

**2534 - 8**  
Winged

**2535 - W8**  
Wingless

**For Square Lower Molars**

**For Round Upper Molars**



**2536 - 8A**  
Winged

**2537 - W8A**  
Wingless

**2538 - 12A**  
Winged

**2539 - 13A**  
Wingless

**For Molar Roots**

**For Partially Broken Molars**



**2540 - 14**  
Winged

**2541 - 14T**  
Edged Winged

**2542 - W14**  
Wingless

**2543 - 14A**  
Winged

**2544 - W14A**  
Wingless

**For Round Molars, Partially Broken, Undersized or Inclined**

**2545 - 26**  
Wingless

**2246 - 51**  
Wingless

**2547 - 202**  
Winged

**2548 - 205**  
Winged